

Antibodies Product Datasheet

Product Name:	MYBPC1 Antibody
Catalog Number:	C16778
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYBPC1. AA range:218-267
Concentration:	1 mg/ml
Gene Name:	MYBPC1
Protein Name:	Myosin-binding protein C slow-type
Other Name:	MYBPC1; MYBPCS; Myosin-binding protein C; slow-type; Slow MyBP-C; C-protein, skeletal muscle slow isoform
Human Swiss Prot No.:	Q00872
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MYBPC1 Polyclonal Antibody detects endogenous levels of MYBPC1 protein.
Storage and Stability:	-20°C/1 year