

Antibodies Product Datasheet

Product Name:	Myosin IXb Antibody
Catalog Number:	C16787
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MYO9B. AA range:304-353
Concentration:	1 mg/ml
Gene Name:	MYO9B
Protein Name:	Unconventional myosin-IXb
Other Name:	MYO9B; MYR5; Unconventional myosin-IXb; Unconventional myosin-9b
Human Swiss Prot No.:	Q13459
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Myosin IXb Polyclonal Antibody detects endogenous levels of Myosin IXb protein.
Storage and Stability:	-20°C/1 year