

Antibodies Product Datasheet

Product Name:	NHERF-2 Antibody
Catalog Number:	C16800
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from NHERF-2 . at AA range: 20-100
Concentration:	1 mg/ml
Gene Name:	SLC9A3R2
Protein Name:	Na(+)/H(+) exchange regulatory cofactor NHE-RF2
Other Name:	SLC9A3R2; NHERF2; Na(+)/H(+) exchange regulatory cofactor NHE-RF2; NHERF-2; NHE3 kinase A regulatory protein E3KARP; SRY-interacting protein 1; SIP-1; Sodium-hydrogen exchanger regulatory factor 2; Solute carrier family 9 isoform A3 regulatory factor 2; T
Human Swiss Prot No.:	Q15599
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NHERF-2 Polyclonal Antibody detects endogenous levels of NHERF-2 protein.
Storage and Stability:	-20°C/1 year