

Antibodies Product Datasheet

Product Name:	NDUFA9 Antibody
Catalog Number:	C16827
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NDUFA9. AA range:87-136
Concentration:	1 mg/ml
Gene Name:	NDUFA9
Protein Name:	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 9 mitochondrial
Other Name:	NDUFA9; NDUF52L; NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 9; mitochondrial; Complex I-39kD; CI-39kD; NADH-ubiquinone oxidoreductase 39 kDa subunit
Human Swiss Prot No.:	Q16795
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NDUFA9 Polyclonal Antibody detects endogenous levels of NDUFA9 protein.
Storage and Stability:	-20°C/1 year