

Antibodies Product Datasheet

Product Name:	NDUFS6 Antibody
Catalog Number:	C16848
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NDUFS6. AA range:75-124
Concentration:	1 mg/ml
Gene Name:	NDUFS6
Protein Name:	NADH dehydrogenase [ubiquinone] iron-sulfur protein 6 mitochondrial
Other Name:	NDUFS6; NADH dehydrogenase [ubiquinone] iron-sulfur protein 6; mitochondrial; Complex I-13kD-A; CI-13kD-A; NADH-ubiquinone oxidoreductase 13 kDa-A subunit
Human Swiss Prot No.:	O75380
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NDUFS6 Polyclonal Antibody detects endogenous levels of NDUFS6 protein.
Storage and Stability:	-20°C/1 year