

Antibodies Product Datasheet

Product Name:	LTB4DH Antibody
Catalog Number:	C16862
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PTGR1. AA range:75-124
Concentration:	1 mg/ml
Gene Name:	PTGR1
Protein Name:	Prostaglandin reductase 1
Other Name:	PTGR1; LTB4DH; Prostaglandin reductase 1; PRG-1; 15-oxoprostaglandin 13-reductase; NADP-dependent leukotriene B4 12-hydroxydehydrogenase
Human Swiss Prot No.:	Q14914
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LTB4DH Polyclonal Antibody detects endogenous levels of LTB4DH protein.
Storage and Stability:	-20°C/1 year