

Antibodies Product Datasheet

Product Name:	Mox1 Antibody
Catalog Number:	C16868
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NOX1. AA range:436-485
Concentration:	1 mg/ml
Gene Name:	NOX1
Protein Name:	NADPH oxidase 1
Other Name:	NOX1; MOX1; NOH1; NADPH oxidase 1; NOX-1; Mitogenic oxidase 1; MOX-1; NADH/NADPH mitogenic oxidase subunit P65-MOX; NOH-1
Human Swiss Prot No.:	Q9Y5S8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mox1 Polyclonal Antibody detects endogenous levels of Mox1 protein.
Storage and Stability:	-20°C/1 year