

Antibodies Product Datasheet

Product Name:	NBEAL1 Antibody
Catalog Number:	C16920
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NBEAL1. AA range:566-615
Concentration:	1 mg/ml
Gene Name:	NBEAL1
Protein Name:	Neurobeachin-like protein 1
Other Name:	NBEAL1; ALS2CR16; ALS2CR17; Neurobeachin-like protein 1; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 16 protein; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 17 protein
Human Swiss Prot No.:	Q6ZS30
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NBEAL1 Polyclonal Antibody detects endogenous levels of NBEAL1 protein.
Storage and Stability:	-20°C/1 year