

Antibodies Product Datasheet

Product Name:	Neuron Navigator 1 Antibody
Catalog Number:	C16979
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Neuron Navigator 1 . at AA range: 510-590
Concentration:	1 mg/ml
Gene Name:	NAV1
Protein Name:	Neuron navigator 1
Other Name:	NAV1; KIAA1151; KIAA1213; POMFIL3; STEERIN1; Neuron navigator 1; Pore membrane and/or filament-interacting-like protein 3; Steerin-1; Unc-53 homolog 1; unc53H1
Human Swiss Prot No.:	Q8NEY1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Neuron Navigator 1 Polyclonal Antibody detects endogenous levels of Neuron Navigator 1 protein.
Storage and Stability:	-20°C/1 year