

Antibodies Product Datasheet

Product Name:	Neuron Navigator 2 Antibody
Catalog Number:	C16980
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Neuron Navigator 2 . at AA range: 2040-2120
Concentration:	1 mg/ml
Gene Name:	NAV2
Protein Name:	Neuron navigator 2
Other Name:	NAV2; HELAD1; KIAA1419; POMFIL2; RAINB1; STEERIN2; Neuron navigator 2; Helicase APC down-regulated 1; Pore membrane and/or filament-interacting-like protein 2; Retinoic acid inducible in neuroblastoma 1; Steerin-2; Unc-53 homolog 2; unc53H2
Human Swiss Prot No.:	Q8IVL1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Neuron Navigator 2 Polyclonal Antibody detects endogenous levels of Neuron Navigator 2 protein.
Storage and Stability:	-20°C/1 year

C16980 | Neuron Navigator 2 Antibody

Tech@assaybiotech.com | Order@assaybiotech.com

Research use only not intended for therapeutic/diagnostic procedures

Last updated December 2022