

Antibodies Product Datasheet

Product Name:	M6A Antibody
Catalog Number:	C16995
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from M6A . at AA range: 140-220
Concentration:	1 mg/ml
Gene Name:	GPM6A
Protein Name:	Neuronal membrane glycoprotein M6-a
Other Name:	GPM6A; M6A; Neuronal membrane glycoprotein M6-a; M6a
Human Swiss Prot No.:	P51674
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	M6A Polyclonal Antibody detects endogenous levels of M6A protein.
Storage and Stability:	-20°C/1 year