

Antibodies Product Datasheet

Product Name:	Glucosidase II α Antibody
Catalog Number:	C17041
Reactivity:	Human, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GANAB. AA range:242-291
Concentration:	1 mg/ml
Gene Name:	GANAB
Protein Name:	Neutral alpha-glucosidase AB
Other Name:	GANAB; G2AN; KIAA0088; Neutral alpha-glucosidase AB; Alpha-glucosidase 2; Glucosidase II subunit alpha
Human Swiss Prot No.:	Q14697
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Glucosidase II α Polyclonal Antibody detects endogenous levels of Glucosidase II α protein.
Storage and Stability:	-20°C/1 year