

Antibodies Product Datasheet

Product Name:	DFNA5 Antibody
Catalog Number:	C17063
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from DFNA5 . at AA range: 200-280
Concentration:	1 mg/ml
Gene Name:	DFNA5
Protein Name:	Non-syndromic hearing impairment protein 5
Other Name:	DFNA5; ICERE1; Non-syndromic hearing impairment protein 5; Inversely correlated with estrogen receptor expression 1; ICERE-1
Human Swiss Prot No.:	O60443
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DFNA5 Polyclonal Antibody detects endogenous levels of DFNA5 protein.
Storage and Stability:	-20°C/1 year