

Antibodies Product Datasheet

Product Name:	NARF Antibody
Catalog Number:	C17103
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NARF. AA range:262-311
Concentration:	1 mg/ml
Gene Name:	NARF
Protein Name:	Nuclear prelamin A recognition factor
Other Name:	NARF; Nuclear prelamin A recognition factor; Iron-only hydrogenase-like protein 2; IOP2
Human Swiss Prot No.:	Q9UHQ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NARF Polyclonal Antibody detects endogenous levels of NARF protein.
Storage and Stability:	-20°C/1 year