

Antibodies Product Datasheet

Product Name:	CIP29 Antibody
Catalog Number:	C17106
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human HCC1. AA range:147-196
Concentration:	1 mg/ml
Gene Name:	SARNP
Protein Name:	SAP domain-containing ribonucleoprotein
Other Name:	SARNP; HCC1; HSPC316; SAP domain-containing ribonucleoprotein; Cytokine-induced protein of 29 kDa; Nuclear protein Hcc-1; Proliferation-associated cytokine-inducible protein CIP29
Human Swiss Prot No.:	P82979
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CIP29 Polyclonal Antibody detects endogenous levels of CIP29 protein.
Storage and Stability:	-20°C/1 year