

Antibodies Product Datasheet

Product Name:	ARA70 Antibody
Catalog Number:	C17114
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ARA70 . at AA range: 320-400
Concentration:	1 mg/ml
Gene Name:	NCOA4
Protein Name:	Nuclear receptor coactivator 4
Other Name:	NCOA4; ARA70; ELE1; RFG; Nuclear receptor coactivator 4; NCoA-4; Androgen receptor coactivator 70 kDa protein; 70 kDa AR-activator; 70 kDa androgen receptor coactivator; Androgen receptor-associated protein of 70 kDa; Ret-activating protein
Human Swiss Prot No.:	Q13772
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ARA70 Polyclonal Antibody detects endogenous levels of ARA70 protein.
Storage and Stability:	-20°C/1 year