

Antibodies Product Datasheet

Product Name:	NCoA-5 Antibody
Catalog Number:	C17115
Reactivity:	Human, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NCOA5. AA range:345-394
Concentration:	1 mg/ml
Gene Name:	NCOA5
Protein Name:	Nuclear receptor coactivator 5
Other Name:	NCOA5; KIAA1637; Nuclear receptor coactivator 5; NCoA-5; Coactivator independent of AF-2; CIA
Human Swiss Prot No.:	Q9HCD5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NCoA-5 Polyclonal Antibody detects endogenous levels of NCoA-5 protein.
Storage and Stability:	-20°C/1 year