

Antibodies Product Datasheet

Product Name:	NOC3L Antibody
Catalog Number:	C17147
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NOC3L. AA range:415-464
Concentration:	1 mg/ml
Gene Name:	NOC3L
Protein Name:	Nucleolar complex protein 3 homolog
Other Name:	NOC3L; AD24; C10orf117; FAD24; Nucleolar complex protein 3 homolog; NOC3 protein homolog; Factor for adipocyte differentiation 24; NOC3-like protein; Nucleolar complex-associated protein 3-like protein
Human Swiss Prot No.:	Q8WTT2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NOC3L Polyclonal Antibody detects endogenous levels of NOC3L protein.
Storage and Stability:	-20°C/1 year