

Antibodies Product Datasheet

Product Name:	PEPT1 Antibody
Catalog Number:	C17548
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SLC15A1. AA range:225-274
Concentration:	1 mg/ml
Gene Name:	SLC15A1
Protein Name:	Solute carrier family 15 member 1
Other Name:	SLC15A1; PEPT1; Solute carrier family 15 member 1; Intestinal H(+)/peptide cotransporter; Oligopeptide transporter; small intestine isoform; Peptide transporter 1
Human Swiss Prot No.:	P46059
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PEPT1 Polyclonal Antibody detects endogenous levels of PEPT1 protein.
Storage and Stability:	-20°C/1 year