

Antibodies Product Datasheet

Product Name:	DHHC-15 Antibody
Catalog Number:	C17595
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZDHHC15. AA range:288-337
Concentration:	1 mg/ml
Gene Name:	ZDHHC15
Protein Name:	Palmitoyltransferase ZDHHC15
Other Name:	ZDHHC15; Palmitoyltransferase ZDHHC15; Zinc finger DHHC domain-containing protein 15; DHHC-15
Human Swiss Prot No.:	Q96MV8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DHHC-15 Polyclonal Antibody detects endogenous levels of DHHC-15 protein.
Storage and Stability:	-20°C/1 year