

Antibodies Product Datasheet

Product Name:	DHHC-9 Antibody
Catalog Number:	C17600
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ZDHHC9. AA range:315-364
Concentration:	1 mg/ml
Gene Name:	ZDHHC9
Protein Name:	Palmitoyltransferase ZDHHC9
Other Name:	ZDHHC9; CXorf11; ZDHHC10; ZNF379; ZNF380; CGI-89; Palmitoyltransferase ZDHHC9; Zinc finger DHHC domain-containing protein 9; DHHC-9; DHHC9; Zinc finger protein 379; Zinc finger protein 380
Human Swiss Prot No.:	Q9Y397
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DHHC-9 Polyclonal Antibody detects endogenous levels of DHHC-9 protein.
Storage and Stability:	-20°C/1 year