

Antibodies Product Datasheet

Product Name:	DECR2 Antibody
Catalog Number:	C17620
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DECR2. AA range:217-266
Concentration:	1 mg/ml
Gene Name:	DECR2
Protein Name:	Peroxisomal 2,4-dienoyl-CoA reductase
Other Name:	DECR2; PDCR; Peroxisomal 2; 4-dienoyl-CoA reductase; pDCR; 2,4-dienoyl-CoA reductase 2
Human Swiss Prot No.:	Q9NUI1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DECR2 Polyclonal Antibody detects endogenous levels of DECR2 protein.
Storage and Stability:	-20°C/1 year