

Antibodies Product Datasheet

Product Name:	LONP2 Antibody
Catalog Number:	C17628
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LONP2. AA range:780-829
Concentration:	1 mg/ml
Gene Name:	LONP2
Protein Name:	Lon protease homolog 2 peroxisomal
Other Name:	LONP2; LONP; Lon protease homolog 2; peroxisomal; Lon protease-like protein 2; Lon protease 2; Peroxisomal Lon protease
Human Swiss Prot No.:	Q86WA8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LONP2 Polyclonal Antibody detects endogenous levels of LONP2 protein.
Storage and Stability:	-20°C/1 year