

Antibodies Product Datasheet

Product Name:	Peroxin 11 β Antibody
Catalog Number:	C17630
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PEX11B. AA range:91-140
Concentration:	1 mg/ml
Gene Name:	PEX11B
Protein Name:	Peroxisomal membrane protein 11B
Other Name:	PEX11B; Peroxisomal membrane protein 11B; Peroxin-11B; Peroxisomal biogenesis factor 11B; Protein PEX11 homolog beta; PEX11-beta
Human Swiss Prot No.:	O96011
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Peroxin 11 β Polyclonal Antibody detects endogenous levels of Peroxin 11 β protein.
Storage and Stability:	-20°C/1 year