

Antibodies Product Datasheet

Product Name:	INPP5J Antibody
Catalog Number:	C17667
Reactivity:	Human, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from INPP5J . at AA range: 850-930
Concentration:	1 mg/ml
Gene Name:	INPP5J
Protein Name:	Phosphatidylinositol 4,5-bisphosphate 5-phosphatase A
Other Name:	INPP5J; PIB5PA; PIPP; Phosphatidylinositol 4; 5-bisphosphate 5-phosphatase A; Inositol polyphosphate 5-phosphatase J
Human Swiss Prot No.:	Q15735
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	INPP5J Polyclonal Antibody detects endogenous levels of INPP5J protein.
Storage and Stability:	-20°C/1 year