

Antibodies Product Datasheet

Product Name:	PC-PLD4 Antibody
Catalog Number:	C17717
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLD4. AA range:457-506
Concentration:	1 mg/ml
Gene Name:	PLD4
Protein Name:	Phospholipase D4
Other Name:	PLD4; C14orf175; Phospholipase D4; PLD 4; Choline phosphatase 4; Phosphatidylcholine-hydrolyzing phospholipase D4
Human Swiss Prot No.:	Q96BZ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PC-PLD4 Polyclonal Antibody detects endogenous levels of PC-PLD4 protein.
Storage and Stability:	-20°C/1 year