

Antibodies Product Datasheet

Product Name:	PFAAP5 Antibody
Catalog Number:	C17725
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human N4BP2L2. AA range:358-407
Concentration:	1 mg/ml
Gene Name:	N4BP2L2
Protein Name:	NEDD4-binding protein 2-like 2
Other Name:	N4BP2L2; CG005; PFAAP5; NEDD4-binding protein 2-like 2; Phosphonoformate immuno-associated protein 5
Human Swiss Prot No.:	Q92802
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PFAAP5 Polyclonal Antibody detects endogenous levels of PFAAP5 protein.
Storage and Stability:	-20°C/1 year