

Antibodies Product Datasheet

Product Name:	C1INH Antibody
Catalog Number:	C17738
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human SERPING1. AA range:342-391
Concentration:	1 mg/ml
Gene Name:	SERPING1
Protein Name:	Plasma protease C1 inhibitor
Other Name:	SERPING1; C1IN; C1NH; Plasma protease C1 inhibitor; C1 Inh; C1Inh; C1 esterase inhibitor; C1-inhibiting factor; Serpin G1
Human Swiss Prot No.:	P05155
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C1INH Polyclonal Antibody detects endogenous levels of C1INH protein.
Storage and Stability:	-20°C/1 year