

Antibodies Product Datasheet

Product Name:	PABP3 Antibody
Catalog Number:	C17765
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000.IHC-p:1:50-300 ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PABPC3. AA range:352-401
Concentration:	1 mg/ml
Gene Name:	PABPC3
Protein Name:	Polyadenylate-binding protein 3
Other Name:	PABPC3; PABP3; PABPL3; Polyadenylate-binding protein 3; PABP-3; Poly(A)-binding protein 3; Testis-specific poly(A)-binding protein
Human Swiss Prot No.:	Q9H361
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PABP3 Polyclonal Antibody detects endogenous levels of PABP3 protein.
Storage and Stability:	-20°C/1 year