

Antibodies Product Datasheet

Product Name:	ATP1AL1 Antibody
Catalog Number:	C17835
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ATP1AL1 . at AA range: 380-460
Concentration:	1 mg/ml
Gene Name:	ATP12A
Protein Name:	Potassium-transporting ATPase alpha chain 2
Other Name:	ATP12A; ATP1AL1; Potassium-transporting ATPase alpha chain 2; Non-gastric H(+)/K(+) ATPase subunit alpha; Proton pump
Human Swiss Prot No.:	P54707
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ATP1AL1 Polyclonal Antibody detects endogenous levels of ATP1AL1 protein.
Storage and Stability:	-20°C/1 year