

Antibodies Product Datasheet

Product Name:	Cdc40 Antibody
Catalog Number:	C17847
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CDC40. AA range:179-228
Concentration:	1 mg/ml
Gene Name:	CDC40
Protein Name:	Pre-mRNA-processing factor 17
Other Name:	CDC40; EHB3; PRP17; PRPF17; Pre-mRNA-processing factor 17; Cell division cycle 40 homolog; EH-binding protein 3; Ehb3; PRP17 homolog; hPRP17
Human Swiss Prot No.:	O60508
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cdc40 Polyclonal Antibody detects endogenous levels of Cdc40 protein.
Storage and Stability:	-20°C/1 year