

Antibodies Product Datasheet

Product Name:	PRP6 Antibody
Catalog Number:	C17850
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), WB, ELISA
Dilutions:	WB 1:500-2000 Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PRPF6. AA range:747-796
Concentration:	1 mg/ml
Gene Name:	PRPF6
Protein Name:	Pre-mRNA-processing factor 6
Other Name:	PRPF6; C20orf14; Pre-mRNA-processing factor 6; Androgen receptor N-terminal domain-transactivating protein 1; ANT-1; PRP6 homolog; U5 snRNP-associated 102 kDa protein; U5-102 kDa protein
Human Swiss Prot No.:	O94906
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PRP6 Polyclonal Antibody detects endogenous levels of PRP6 protein.
Storage and Stability:	-20°C/1 year