

Antibodies Product Datasheet

Product Name:	PRPF18 Antibody
Catalog Number:	C17852
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PRPF18. AA range:182-231
Concentration:	1 mg/ml
Gene Name:	PRPF18
Protein Name:	Pre-mRNA-splicing factor 18
Other Name:	PRPF18; HPRP18; Pre-mRNA-splicing factor 18; PRP18 homolog; hPRP18
Human Swiss Prot No.:	Q99633
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PRPF18 Polyclonal Antibody detects endogenous levels of PRPF18 protein.
Storage and Stability:	-20°C/1 year