

Antibodies Product Datasheet

Product Name:	DDX52 Antibody
Catalog Number:	C17865
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DDX52. AA range:550-599
Concentration:	1 mg/ml
Gene Name:	DDX52
Protein Name:	Probable ATP-dependent RNA helicase DDX52
Other Name:	DDX52; ROK1; HUSSY-19; Probable ATP-dependent RNA helicase DDX52; ATP-dependent RNA helicase ROK1-like; DEAD box protein 52
Human Swiss Prot No.:	Q9Y2R4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	DDX52 Polyclonal Antibody detects endogenous levels of DDX52 protein.
Storage and Stability:	-20°C/1 year