

Antibodies Product Datasheet

Product Name:	FAM35A/B Antibody
Catalog Number:	C17989
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FAM35A/B. AA range:394-443
Concentration:	1 mg/ml
Gene Name:	FAM35A/FAM35B
Protein Name:	Protein FAM35A/B
Other Name:	FAM35A; Protein FAM35A; FAM35B; Protein FAM35B
Human Swiss Prot No.:	Q86V20/A6NDK8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FAM35A/B Polyclonal Antibody detects endogenous levels of FAM35A/B protein.
Storage and Stability:	-20°C/1 year