

Antibodies Product Datasheet

Product Name:	FT β Antibody
Catalog Number:	C17991
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FNTB. AA range:255-304
Concentration:	1 mg/ml
Gene Name:	FNTB
Protein Name:	Protein farnesyltransferase subunit beta
Other Name:	FNTB; Protein farnesyltransferase subunit beta; FTase-beta; CAAX farnesyltransferase subunit beta; Ras proteins prenyltransferase subunit beta
Human Swiss Prot No.:	P49356
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FT β Polyclonal Antibody detects endogenous levels of FT β protein.
Storage and Stability:	-20°C/1 year