

Antibodies Product Datasheet

Product Name:	TIMAP Antibody
Catalog Number:	C18003
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP1R16B. AA range:364-413
Concentration:	1 mg/ml
Gene Name:	PPP1R16B
Protein Name:	Protein phosphatase 1 regulatory inhibitor subunit 16B
Other Name:	PPP1R16B; ANKRD4; KIAA0823; Protein phosphatase 1 regulatory inhibitor subunit 16B; Ankyrin repeat domain-containing protein 4; CAAX box protein TIMAP; TGF-beta-inhibited membrane-associated protein; hTIMAP
Human Swiss Prot No.:	Q96T49
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TIMAP Polyclonal Antibody detects endogenous levels of TIMAP protein.
Storage and Stability:	-20°C/1 year