

Antibodies Product Datasheet

Product Name:	PPP1R3A Antibody
Catalog Number:	C18014
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPP1R3A. AA range:647-696
Concentration:	1 mg/ml
Gene Name:	PPP1R3A
Protein Name:	Protein phosphatase 1 regulatory subunit 3A
Other Name:	PPP1R3A; PP1G; Protein phosphatase 1 regulatory subunit 3A; Protein phosphatase 1 glycogen-associated regulatory subunit; Protein phosphatase type-1 glycogen targeting subunit; RG1
Human Swiss Prot No.:	Q16821
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PPP1R3A Polyclonal Antibody detects endogenous levels of PPP1R3A protein.
Storage and Stability:	-20°C/1 year