

Antibodies Product Datasheet

Product Name:	PP2Cε Antibody
Catalog Number:	C18028
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PPM1L. AA range:119-168
Concentration:	1 mg/ml
Gene Name:	PPM1L
Protein Name:	Protein phosphatase 1L
Other Name:	PPM1L; PP2CE; Protein phosphatase 1L; Protein phosphatase 1-like; Protein phosphatase 2C isoform epsilon; PP2C-epsilon
Human Swiss Prot No.:	Q5SGD2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PP2Cε Polyclonal Antibody detects endogenous levels of PP2Cε protein.
Storage and Stability:	-20°C/1 year