

Antibodies Product Datasheet

Product Name:	Calgizzarin Antibody
Catalog Number:	C18040
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Calgizzarin . at AA range: 30-110
Concentration:	1 mg/ml
Gene Name:	S100A11
Protein Name:	Protein S100-A11
Other Name:	S100A11; MLN70; S100C; Protein S100-A11; Calgizzarin; Metastatic lymph node gene 70 protein; MLN 70; Protein S100-C; S100 calcium-binding protein A11
Human Swiss Prot No.:	P31949
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Calgizzarin Polyclonal Antibody detects endogenous levels of Calgizzarin protein.
Storage and Stability:	-20°C/1 year