

Antibodies Product Datasheet

Product Name:	APOBEC4 Antibody
Catalog Number:	C18080
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human APOBEC4. AA range:309-358
Concentration:	1 mg/ml
Gene Name:	APOBEC4
Protein Name:	Putative C->U-editing enzyme APOBEC-4
Other Name:	APOBEC4; C1orf169; Putative C->U-editing enzyme APOBEC-4; Apolipoprotein B mRNA-editing enzyme catalytic polypeptide-like 4
Human Swiss Prot No.:	Q8WW27
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	APOBEC4 Polyclonal Antibody detects endogenous levels of APOBEC4 protein.
Storage and Stability:	-20°C/1 year