

Antibodies Product Datasheet

Product Name:	Ran BP-6 Antibody
Catalog Number:	C18159
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RANBP6. AA range:800-849
Concentration:	1 mg/ml
Gene Name:	RANBP6
Protein Name:	Ran-binding protein 6
Other Name:	RANBP6; Ran-binding protein 6; RanBP6
Human Swiss Prot No.:	O60518
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ran BP-6 Polyclonal Antibody detects endogenous levels of Ran BP-6 protein.
Storage and Stability:	-20°C/1 year