

Antibodies Product Datasheet

Product Name:	RGPD2/3/4/6/8 Antibody
Catalog Number:	C18163
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RGPD2/3/4/6/8. AA range:577-626
Concentration:	1 mg/ml
Gene Name:	RGPD8
Protein Name:	RANBP2-like and GRIP domain-containing protein 8
Other Name:	
Human Swiss Prot No.:	O14715/Q53T03/Q7Z3J3/A6NKT7/P0C839
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	RGPD2/3/4/6/8 Polyclonal Antibody detects endogenous levels of RGPD2/3/4/6/8 protein.
Storage and Stability:	-20°C/1 year