

Antibodies Product Datasheet

Product Name:	Epac Antibody
Catalog Number:	C18168
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAPGEF3. AA range:712-761
Concentration:	1 mg/ml
Gene Name:	RAPGEF3
Protein Name:	Rap guanine nucleotide exchange factor 3
Other Name:	RAPGEF3; CGEF1; EPAC; EPAC1; Rap guanine nucleotide exchange factor 3; Exchange factor directly activated by cAMP 1; Exchange protein directly activated by cAMP 1; EPAC 1; Rap1 guanine-nucleotide-exchange factor directly activated by cAMP;
Human Swiss Prot No.:	O95398
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Epac Polyclonal Antibody detects endogenous levels of Epac protein.
Storage and Stability:	-20°C/1 year