

Antibodies Product Datasheet

Product Name:	GAP1-InsP4 BP Antibody
Catalog Number:	C18188
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RASA3. AA range:247-296
Concentration:	1 mg/ml
Gene Name:	RASA3
Protein Name:	Ras GTPase-activating protein 3
Other Name:	RASA3; Ras GTPase-activating protein 3; GAP1(IP4BP); Ins P4-binding protein
Human Swiss Prot No.:	Q14644
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GAP1-InsP4 BP Polyclonal Antibody detects endogenous levels of GAP1-InsP4 BP protein.
Storage and Stability:	-20°C/1 year