

Antibodies Product Datasheet

Product Name:	G3BP2 Antibody
Catalog Number:	C18193
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human G3BP2. AA range:433-482
Concentration:	1 mg/ml
Gene Name:	G3BP2
Protein Name:	Ras GTPase-activating protein-binding protein 2
Other Name:	G3BP2; KIAA0660; Ras GTPase-activating protein-binding protein 2; G3BP-2; GAP SH3 domain-binding protein 2
Human Swiss Prot No.:	Q9UN86
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	G3BP2 Polyclonal Antibody detects endogenous levels of G3BP2 protein.
Storage and Stability:	-20°C/1 year