

Antibodies Product Datasheet

Product Name:	IQGAP1 Antibody
Catalog Number:	C18194
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human IQGAP1. AA range:247-296
Concentration:	1 mg/ml
Gene Name:	IQGAP1
Protein Name:	Ras GTPase-activating-like protein IQGAP1
Other Name:	IQGAP1; KIAA0051; Ras GTPase-activating-like protein IQGAP1; p195
Human Swiss Prot No.:	P46940
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IQGAP1 Polyclonal Antibody detects endogenous levels of IQGAP1 protein.
Storage and Stability:	-20°C/1 year