

Antibodies Product Datasheet

Product Name:	Rab 35 Antibody
Catalog Number:	C18245
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human RAB35. AA range:99-148
Concentration:	1 mg/ml
Gene Name:	RAB35
Protein Name:	Ras-related protein Rab-35
Other Name:	RAB35; RAB1C; RAY; Ras-related protein Rab-35; GTP-binding protein RAY; Ras-related protein Rab-1C
Human Swiss Prot No.:	Q15286
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Rab 35 Polyclonal Antibody detects endogenous levels of Rab 35 protein.
Storage and Stability:	-20°C/1 year