

Antibodies Product Datasheet

Product Name:	PTP IA-2 β Antibody
Catalog Number:	C18297
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PTPRN2. AA range:206-255
Concentration:	1 mg/ml
Gene Name:	PTPRN2
Protein Name:	Receptor-type tyrosine-protein phosphatase N2
Other Name:	PTPRN2; KIAA0387; Receptor-type tyrosine-protein phosphatase N2; R-PTP-N2; Islet cell autoantigen-related protein; IAR; ICAAR; Phogrin
Human Swiss Prot No.:	Q92932
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PTP IA-2 β Polyclonal Antibody detects endogenous levels of PTP IA-2 β protein.
Storage and Stability:	-20°C/1 year