

Antibodies Product Datasheet

Product Name:	PTP ζ Antibody
Catalog Number:	C18303
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PTPRZ1. AA range:122-171
Concentration:	1 mg/ml
Gene Name:	PTPRZ1
Protein Name:	Receptor-type tyrosine-protein phosphatase zeta
Other Name:	PTPRZ1; HTPZP2; PTPRZ; PTPRZ2; PTPZ; Receptor-type tyrosine-protein phosphatase zeta; R-PTP-zeta; Protein-tyrosine phosphatase receptor type Z polypeptide 1; Protein-tyrosine phosphatase receptor type Z polypeptide 2; R-PTP-zeta-2
Human Swiss Prot No.:	P23471
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PTP ζ Polyclonal Antibody detects endogenous levels of PTP ζ protein.
Storage and Stability:	-20°C/1 year